


VISCO JET® CRACK - 120 mm

P/N: 509-17120-00



Particularly in the case of applications with very viscous media such as pastes and gels, a very reliable stirring process in the entire vessel without air entry is important. This is what the innovative stirring system VISCO JET® offers!

The particular shape of the blades triggers an intrinsically dynamic flow, which shortens the process times with significantly better mixing results compared to conventional impellers. An all-rounder for aqueous and highly viscous substances. You can find out more about this [here](#).

The VISCO JET® Crack 120 is ideal for vessels with the following characteristics:

- ideal vessel diameter: 200- 300 mm
- ideal ground clearance: 36 mm
- ideal filling level: 200- 300 mm
- ideal volume: 6,300 - 21,200 ml

The VISCO JET® Crack 120 is made of stainless steel

A stirrer shaft with a diameter of 10 mm is included in the scope of delivery.

VISCO JET® CRACK - 120 mm - Teknik veri

Blade ölçüsü	120 mm
Materyal	stainless steel V4A
Uzunluk	500 mm
Şaft çapı	10 mm
Vesel çapı için	170 - 300 mm
Rotasyon hız derecesi	120 - 500 rpm