

Hei-FLOW Core 600

P/N: 523-50060-00



This pump of the Hei-Flow Core series is particularly suitable when high volumes need to be transported easily. A flow rate of 4.0 - 4,151 ml / min can be set for this purpose. This model is designed for the use of single-channel pump heads.

Even with a load change, the speed of 50 – 600 rpm is kept constant at a control accuracy of $\pm 2\%$. It is easy to switch to the reverse direction of rotation.

Hei-FLOW Core 600 - Teknik veri

İzin verilen ortam koşulları	5 - 31 °C at 80% rel. humidity 32 - 40 °C decreasing linearly up to max. 50% rel. humidity
Ağırlık	7,1 kg
Protection class (EN 60529)	IP 54
Güç sağlamak	100 W
Dimensions (w/d/h)	166 x 256 x 225 mm
Akış derecesi tekli kanal pompaları	4.0 – 4,151 ml/min
Akış derecesi- çoklu kanal pompaları	-
Akış derece göstergesi	±5 %
Hız derecesi	50 - 600 rpm
Hız ayarları	scale
Elektronik hız kontrol	digital
Kontrol doğruluğu motor	±2 %
Dönme yönünü seçin	CW / CCW
Motor gücü	100 W
Analog arayüz	-
Dijital arayüz	-
Akış derece göstergesi	-
Hacim dozlama	-
Entegre vakum kontrol	-
Düzgün başlama	-
Elektronik fren	-
Ayak-pedal bağlantısı	-
Sürekli operasyon(24/7)	yes
Güvenli özellik	overheat protection