



## Hei-VAP Core - 3000

Display	3.5" LCD
Humidity Control	5-31% RH (5-80%) 32-40% RH (50-55%)
Compliance	DIN EN 60529
Protection	IP 20
Material	Stainless Steel
Control	Touch
Height	155 mm
Power	1,300 W
Temperature Range	20 – 210 °C
Temperature Accuracy	±1K
Temperature Sensor	PT 1000 Class B, 5% Accuracy
Display	3.5" LCD
Control	Touch
Temperature	250 °C
Material	AISI 316L
Height	255 mm
Capacity	4.5 l
Weight	-
Dimensions	-
Material	-
Power	-
Weight	-
Power	1,400 W
Weight	15 kg
Dimensions G3 (w/d/h)	739 x 487 x 845 mm
Dimensions G3 XL (w/d/h)	739 x 487 x 920 mm

## Premium Laboratory Equipment

Area	1,400 cm <sup>2</sup>
XL Area	2,200 cm <sup>2</sup>
Volume	5 l
Speed	10 - 280 rpm
Hei-CHILL	-
Dynamic AUTOaccurate	-
Solvent library	-